



TRAVEL POWER CHAIR

EXPLORE THE WORLD

Whether you are riding the rails or cruising the seas, you need a travel power wheelchair that is lightweight and portable. Consider a Pride travel power wheelchair for your next adventure.



SCAN THE QR CODE
TO LEARN MORE ABOUT
TRAVEL MOBILITY.

PRIDE
MOBILITY
LIVE YOUR BEST®

1GO CHAIR



FEATURES

- Easily disassembles into 5 lightweight pieces
- Dual swing-away, under seat storage bins
- Number one for portability
- Heaviest piece is only 36 lbs.
- Available in 6 colors

ACCESSORIES

- Armrest mount cane clip
- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear mount bracket
- Saddlebag
- Safety flag
- Weather cover
- XLR USB charger
- Additional battery pack

Ruby Red

Pearl White

Amber Orange

Citrine Yellow

Rose Quartz

Sapphire Blue

MILITARY PATCHES



300 lbs.	Up To 3.8 mph	25.5"	Up To 13.2
Weight Capacity	Maximum Speed	Turning Radius	Miles Per Charge at 200 lbs.2

JAZZY® PASSPORT 1



FEATURES

- Compact and lightweight design makes traveling a breeze
- Folds in just a few simple steps
- Convenient, under seat mesh storage bag and cup holder
- Airline safe lithium-ion battery available3
- Only 59.8 lbs. without batteries
- Available in black

ACCESSORIES

- XLR USB charger

OPTIONAL

- Stealth Simplicity G Cushion, 18"x18"



Black

Folded Size
31" x 16" x 23.5"



MILITARY PATCHES



250 lbs.	Up To 3.6 mph	24.25"	Up To 8.6	Up To 7.9
Weight Capacity	Maximum Speed	Turning Radius	Miles Per Charge at 220 lbs.2	Miles Per Charge at 250 lbs.2

GO CHAIR® AND JAZZY® PASSPORT WARRANTY

5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month limited on batteries



AUTHORIZED PROVIDER
704-821-7777

www.ShopMobility.com

HCPCS Code: K0899

1. Pride Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2. Range varies based on user weight, terrain and battery type.
3. Lithium-ion battery is subject to federal regulations and airline policy. Please consult your airline prior to flight.

PRIDE
MOBILITY
LIVE YOUR BEST®